

How it works

Simply back your wheelchair up to the SPEEDY-LIFT, push the 'IN' button and watch the SPEEDY-LIFT pick up your chair and tuck it into place behind the driver's seat.

Features

- A full cycle takes less than 30 seconds
- The compact design allows you to keep the use of up to 6 seats
- Equipped with an automatic locking system that fastens your wheelchair in place

Specifications

- Maximum capacity: 95 lbs
- Designed for rigid frame manual wheelchairs
- Optional adaptor kit available for folding wheelchairs
- Can be paired with the LINK, the ASENTO or the XL-SEAT





SPEEDY-LIFT

Fast wheelchair lift

Simplified access

The SPEEDY-LIFT is designed to fit adjacent to the sliding door in most minivans on the driver or passenger side. Because of its unparalleled speed, precision and ease, the SPEEDY-LIFT is an ideal solution for someone looking for simplified vehicle access.

Quick setup

The SPEEDY-LIFT is easily installed into the second row seat mounts with no structural modifications to your vehicle. It is also easily removed to return the vehicle to its original condition without affecting resale value.

Adaptive equipment

Our products qualify as adaptive equipment and are eligible for reimbursement under the terms of Mobility programs offered by various vehicle manufacturers.

About us



We take our customers' safety very seriously. We are dedicated to manufacturing quality mobility equipment. Our products undergo a strict certification process.

Certification of our products involves different tests such as: crash testing, cycle testing and weight capacity testing. We are proud to offer a line of products that exceed the industry's standards.



Dealer:



145 Damase-Breton | Saint-Lambert, Qc G0S 2W0

TF 1.866.641.0419 | Fax 418.889.9838 | P 418.889.0419

www.adaptsolutions.com



Due to the diversity of mobility solutions available to you, please contact your local mobility dealer for a personalized evaluation and overview of options tailored to your needs.